

Total Exposures: 76 (4 channels, 19 frames)

Data from: LAS X 3.5.2.18963

Dimensions

Dimension	Logical Size	Physical Length	Start Position	End Position	Pixel Size / Voxel Size
Х	512	46.13 μm	0 μm	46.13 μm	0.09 μm
Υ	512	46.13 μm	0 μm	46.13 μm	0.09 μm
Z	19	5.4 μm	1,187.18 µm	1,192.58 µm	0.3 μm

Channels

LUT	Resolution	Min	Max	STED: DetectorMode / Huygens saturation factor / Wavelength
Blue	8	0	255	/ /
Magenta	8	0	255	/ /
Red	8	0	255	/ /
Gray	8	0	255	/ /

Time Stamps:

Frame (Show All)	Relative Time (s)	Absolute Time (h:m:s.ms)	Date
1	0.000	6:07:58 PM.771	6/11/2019
76	186.955	6:11:05 PM.726	6/11/2019

Confocal Settings

Name	Value
Rotator Angle	0°
Scan Mode	xyz
Scan Direction X	Unidirectional
Scan Speed	400 Hz
Version Number	15
StagePosX	56,073.66 μm
StagePosY	37,948.81 μm

ZPosition	1,192.58 μm
IsSuperZ	0
Magnification	63
ObjectiveName	HC PL APO CS2 63x/1.40 OIL
Immersion	OIL
Numerical Aperture	1.4
RefractionIndex	1.518
Zoom	4
Pinhole	95.5 μm
PinholeAiry	999.44 mAU
EmissionWavelength for PinholeAiry Calculation	580 nm
FrameAverage	1
LineAverage	4
FrameAccumulation	1
Line_Accumulation	1
IsUserSettingNameSet	0
IsRoiScanEnable	0

Filter Wheels / Other Motorized Devices

THE THEORY CHIEF INTERESTEE	
Device Name	Filter Name/Position
Excitation Beam Splitter	TD 488/561/633
Galvo Slider	Galvo X Normal
Notch FW 2	Empty
Polarization FW	NF 488/561/633
Galvo Resonant Pan	Galvo X Pan Center
Target Slider	Target Park
X2 Lens Changer	CS2 UV Optics 1

Lasers

LaserName	OutputPower
Diode 405	On
Argon	On, 0.0000 %
DPSS 561	On
HeNe 633	On

Laser Lines

Laser Line	Intensity				
(405 nm)	Shutter: on, Intensity: 1.5065%				
(633 nm)	Shutter: on, Intensity: 1.0005%				

Detectors

Name	Channel	Туре	Location	Active	Gain	Offset	Gate Start	Gate End	Gate Ref. Wavelength
------	---------	------	----------	--------	------	--------	---------------	-------------	-------------------------

PMT 2	Channel 2	PMT (410nm - 527nm)	Internal	Active	750	0	Time Gating not supported
PMT 4	Channel 4	PMT (638nm - 747nm)	Internal	Active	680	0	Time Gating not supported

Sequential Setting Nr.1

PMT 2	Active, Gain: 0, Offset: 0
PMT 4	Active, Gain: 0, Offset: 0
Laser Line (405 nm)	Intensity: 1.5065%
Laser Line (458 nm)	Intensity: 0.0000%
Laser Line (488 nm)	Intensity: 0.0000%
Laser Line (514 nm)	Intensity: 0.0000%
Laser Line (561 nm)	Intensity: 0.0000%
Laser Line (633 nm)	Intensity: 1.0005%

Sequential Setting Nr.2

HyD 3	Active, Gain: 0, Offset: 0
PMT Trans	Active, Gain: 0, Offset: 0
Laser Line (405 nm)	Intensity: 0.0000%
Laser Line (458 nm)	Intensity: 0.0000%
Laser Line (488 nm)	Intensity: 0.0000%
Laser Line (514 nm)	Intensity: 0.0000%
Laser Line (561 nm)	Intensity: 3.0012%
Laser Line (633 nm)	Intensity: 0.0000%